

Alpha Advanced and Alpha Basic



The Alpha models have forearm supports with brake handles at the end of the pads. The Advanced is a joy-stick style handle which is push and pulled. It can be rotated on the ball joint mount to also be pushed side to side for braking and setting the park.

The Basic handle is the Dolomite looped brake handle mounted vertically. You squeeze to brake and push out to set the park.



Alpha models can also be built with one hand brakes on either side, using either style brake handle.



A large pad can be substituted for the individual troughs. Similar to a 3 sided tray, it permits the user to lean in at any point and reach across to the brake handles in front.



Walker base models

We use the Legacy, Maxi or Symphony as our bases to build these custom rollators. The Legacy and Symphony come in 3 seat heights. These models allow us to configure a walker with forearm pads on a smaller scale. The seat height of the super low models is only 18" and 19". The Legacy also has the option for us to use the tall handles to gain 3" additional handle height.

The Maxi+ has a capacity of 440 lbs., a wider seat and comes in 2 seat heights.



To order a custom Dolomite, provide some details about your client, weight, height, ability or special need. We'll design a system and send you a quote.

We would prefer to install these accessories and ship you a completed walker. But with some instruction these can be added to an existing walker.

Contact us:

Phone: 888-347-4537

fax 724-695-2922

email info@clarkehealthcare.com,



Clarke Health Care Products

Let us customize a Dolomite walker for your client's special needs.

Using the Dolomite Legacy, Symphony or Maxi as our base, we can add accessories or modify handles and brakes to accommodate their ability to use a rollator.



Clarke Health Care Products

7830 Steubenville Pike

Oakdale PA 15071

Phone: 888-347-4537

Email: info@clarkehealthcare.com

www.clarkehealthcare.com

One Hand Brakes



For amputees, stroke survivors, one hand allows control of the braking action from one hand. The two brake

cables will exit the brake handle, with one cable crossing over to the other brake pad. The other handle can be a forearm trough or only a handle to grasp.

Crossbar



Another one hand brake option is the Crossbar, this is a tubular bar that replaces the handles and back strap. It resembles the push bar of a shopping cart with a brake handle in the center. The padded bar provides larger area to grasp and push the rollator. The looped brake handle in the center can be a right or left one-hand brake system.

Forearm Troughs



Forearm pads, or troughs allow the client to lean on their forearms with brake control at the end of the pads. The brakes can be on mounted on both sides or controlled from one side. The pads are lightly cushioned and have some adjustment in angle and position in front of the body.

Note that using the forearm pads, you restrict use of the seat. Depending on position and height, the pads may be in the way to comfortably sit on the seat.

HP Handles



HP handles or hemiplegic handles are mounted on the Symphony model. The one hand braking system uses a padded bar on the non-braking side. This bar can be adjusted in height and angle to the body, so that you can push with your body while steering and braking from the other side.

Slow Down Brakes



If the walker seems to be moving too easily or getting away from the client, add Slow Down brakes. They do not affect the braking action, but only add tension or drag on the wheels. The thumb screw will tighten a nylon wheel onto the back tire. The left and right sides can be individually adjusted. For stroke survivors who may be stronger on one side, tighten that side so they walk in a straight line. For gait training, tighten both to slow the tire rotation. You will still have all the braking action from the brake handles. There is a hidden nut under the thumbscrew to lock in your adjustment, if needed.





Custom Dolomite Rollator Worksheet for quote

Mark For _____

Walker	Choose one	4'10"-5'2" suggested height	5'2"-5'5" suggested height	5'4"-6'2" suggested height
<330 lbs Legacy Blue		<input type="checkbox"/> D12054B SuperLow	<input type="checkbox"/> D12052B Low	<input type="checkbox"/> D12050B Standard
<330 lbs Legacy Red		<input type="checkbox"/> D12054R SuperLow	<input type="checkbox"/> D12052R Low	<input type="checkbox"/> D12050R Standard
<275 lbs. Symphony		<input type="checkbox"/> D12162 SuperLow	<input type="checkbox"/> D12161 Low	<input type="checkbox"/> D12160 Standard
<452 lbs. Maxi+			<input type="checkbox"/> D12122 Maxi+ 550	<input type="checkbox"/> D12120 Maxi+ 650

Arm pads	Choose one setup			
D13005 R or L Forearm trough		<input type="checkbox"/> D13005L one pad left	<input type="checkbox"/> D13005R one pad right	<input type="checkbox"/> D13005R + D13005L
D13018 Full 3-sided padded platform with brakes at front edge		<input type="checkbox"/>		

Brake handle	Choose one setup	One-hand Left	One Hand Right	Braking both hands
Advanced, joystick, braking is push/pull		<input type="checkbox"/> D12515 Left mount	<input type="checkbox"/> D12515 Right mount	<input type="checkbox"/> D13051 x 2
Basic looped brake, braking action is squeeze		<input type="checkbox"/> D13026L left mount	<input type="checkbox"/> D13026R right mount	<input type="checkbox"/> D13026R+ D13026L
If using one hand brake do you require a grip handle on other side		<input type="checkbox"/> D13026 + D70057		

Alternate one hand braking systems, will not work with arm pads above				
Crossbar, front push bar w/ one hand brake mounted center, <i>not available for Maxi+</i>		<input type="checkbox"/> D1536813 left mount	<input type="checkbox"/> D1531576 right mount	
HP Handle, one hand brake w/ angle adjustable bar on other side, <i>only available on Symphony</i>		<input type="checkbox"/> D12611R left brake, right HP	<input type="checkbox"/> D12611L right brake, left HP	

Accessories	
D12533 Tray	<input type="checkbox"/>
D13041 Slow Down Brakes	<input type="checkbox"/>
D12606 Cane holder	<input type="checkbox"/>
D12523 Under seat bag	<input type="checkbox"/>
D12682 Oxygen Tank Holder	<input type="checkbox"/>
P12519 Seat Pad	<input type="checkbox"/>
D12686 Allen Wrench tool for pad adjustment	<input type="checkbox"/>

Ordering company _____

Return quote to _____

Email _____

Fax: _____ Ph: _____

This walker will be built in our shop. It cannot be shipped totally assembled. The handles will need to be inserted and adjusted for height and length. It is not recommended for drop ship delivery unless a technician or visiting nurse can be present to make the adjustments to fit the client. There is no option for a backstrap with this assembly. Note, the part numbers in the chart above do not indicate all components necessary to complete the walker assembly.

info@clarkehealthcare.com

888-347-4537