

# BORIS...

toilet and shower chair with tilt-in-space positioning

...we achieve more mobility!



**Basic model BORIS Size 2**  
Also available with headrest and foot supports, arm supports and insert pads





**Bathing and hygiene aids**

The anatomically molded, deep seat pan made of PE plastic with torso and side guides provides particularly high seat stability and tilt-in-space positioning

- Available in five different sizes
- Robust and easy-care sub-frame made of stainless steel
- Height adjustable seat pan for adaptation to the respective constructional situation
- The seat angle is steplessly adjustable via a gas pressure spring from (-) 5° to 30°
- For support in sitting position, the BORIS toilet and shower chair can be combined with the fixation system FLEXI-AKTIV
- We offer head rests and foot supports, insert and arm support pads and a plug-in spray protection as accessories for individual adaptation
- HMV No. 33.40.04.2008



Plug-in arm support for stabilisation and for keeping the upper body upright



Tilt-in-space seat angle is steplessly adjustable via a gas pressure spring from (-) 5° to 30°

Anatomically molded seat pan with extremely good side guides for stabilization of the upper body. Two-part insert pad for additional sitting comfort



Plug-in spray protection available in two sizes



Foldable and angle adjustable foot support made of high-quality ABS material



Height and angle adjustable head rest for support of head control



Practical toilet arrangement

Dimensions inches	Size 0	Size 1	Size 2	Size 3	Size 4
Seat width at front	12	14	16	18	20
Seat width at rear	10	12	14	16	18
Seat depth	12	14	16	18	18.7
Back height	19	20	21.8	23.7	25.3
Seat angle	(-)5 -30°	(-)5 -30°	(-)5 -30°	(-)5 -30°	(-)5 -30°
height over a toilet	18.3-22.2	17.5-21.5	16.3-20.2	16.3-20.2	16.3-20.2
Wide hygiene opening at front	4	4	4	4	4
Wide hygiene opening at back	4.7	4.7	5.3	5.3	5.3
Total width	21.8	23.8	24.5	27.8	28.8
Total length	25.3	25.3	25.3	27.3	27.3
Max. lower leg length	14.8	14.8	18.7	18.7	18.7
Max. load (lb)	132	176	220	220	330

Frame color: Stainless steel



Clarke Health Care Products

www.clarkehealthcare.com 888-347-4537

