

**AQUATEC®**

# Ergo VIP

# Ergo Dual VIP



**Modular personal care solutions  
for complex needs**

ERGONOMIC DESIGN | SAFETY | INDEPENDENCE



# Modular personal care solutions for those with more complex seating and positioning needs

## Highly customizable


The Aquatec Ergo VIP and Ergo Dual VIP are modular tilt in space shower chair commodes designed to suit a wide range of individuals and their postural needs. Both models benefit from a -5° to 40° seat tilt angle, whilst the Dual VIP has an additional backrest recline of 0° to 35° for those with more complex seating and positioning requirements.

As with all Ergo models, the VIP and Dual VIP can be tailored to each individual by making some simple adjustments or by adding some easy to fit, off the shelf accessories. Furthermore, their ergonomic design ensures ease of use for the carer whilst providing optimum safety and comfort to the user.

See full range of accessories online at [www.clarkehealthcare.com](http://www.clarkehealthcare.com)

### NEW Contoured, ergonomic seat

- ▶ Based on the imprints people leave behind when seated to ensure a more upright and stable seating position.
- ▶ A 5° seat tilt to lift the knees into the natural squat position for toileting.
- ▶ Transfer supports at the front to allow the user to position themselves into a comfortable position for toileting and showering.



### NEW Ergonomic tilt

- ▶ The new pivot and slide mechanism means that during tilt, the chair now moves downwards and forwards which keeps the user's center of gravity close to that of the chair. This makes the tilt operation lighter and more intuitive for the carer, creates a smoother and more comfortable experience for the user and has reduced the overall footprint of the product without affecting its stability.



**NEW**

FEATURES

**-5° to 40°**

**1** **Extremely stable**

- ▶ Enhanced stability for additional safety and reassurance

**2** **Ergonomic tilt**


- ▶ Ergonomic tilt with weight shifting abilities

**3** **Innovative ergonomic seat**

- ▶ Contoured, ergonomic seat plate encourages a more upright seating position for enhanced independence and comfort

**4** **Smaller footprint**

- ▶ Optimized wheel base makes it easier to maneuver in small spaces



## Other features



▲ Stainless steel corrosion free frame.



▲ Extremely stable even at maximum tilt.



▲ Machine washable tension adjustable backrest.



▲ Armrest height, footrest height and seat height can be easily adjusted.



5° - 35°

▲ 5° - 35° Backrest Recline

**Ergo Dual VIP Only**



# Meet the Ergo VIP and Ergo Dual VIP ...

---

Modular tilt in space models to suit a wide range of individuals and their postural needs.



***Ergo VIP***

▲ -5°-40° seat tilt angle

---



***Ergo DualVIP***

▲ -5°-40° seat tilt angle with an additional  
backrest recline of 5° - 35°

---

## Technical Data ▼

For more comprehensive information about this product, including the product's user manual, please see your Clarke Health Care rep or Aquatec dealer.



Seat width ▼



Seat depth ▼



Total width ▼



Total depth ▼



Width between armrests ▼

<b>Ergo VIP</b>	18.7"	18.7"	22"	39"	17.75"
<b>Ergo Dual VIP</b>	18.7"	18.7"	22"	39"	17.75"



Tilt angle seat ▼



Tilt angle back ▼



Weight ▼



Load capacity ▼

<b>Ergo VIP</b>	-5° - 40°	-	51.8 lb.	330 lb.
<b>Ergo Dual VIP</b>	-5° - 40°	5° - 35°	62.4 lb.	330 lb.

## Frame color ▼



White

## Regulatory data ▼



**Clarke Healthcare Products**  
 7830 Steubenville Pike, Oakdale, PA 15071  
[www.clarkehealthcare.com](http://www.clarkehealthcare.com) **888-347-4537**

