



Custom Dolomite Rollator Worksheet for quote, return to info@clarkehealthcare.com or fax: 724-695-2922

Walker	Choose one model, size or color	4'10"-5'2" suggested height	5'2"-5'5" suggested height	5'4"-6'2" suggested height
	Legacy , <330 lbs cap. choose frame color	<input type="checkbox"/> Legacy 450 SuperLow Seat height 19", handle 28-31" <input type="checkbox"/> Blue	<input type="checkbox"/> Legacy 520 Low Seat height 22", handle 31-34" <input type="checkbox"/> Red	<input type="checkbox"/> Legacy 600 Standard Seat height 24", handle 33-37"
	Symphony, <275 lbs cap. Silver frame	<input type="checkbox"/> Sym. 450 SuperLow Seat height 19", handle 25-30"	<input type="checkbox"/> Sym. 520 Low Seat height 22", handle 27.5-33"	<input type="checkbox"/> Sym. 600 Standard Seat height 24", handle 30-36"
	Maxi+, <452 lbs cap., Blue frame		<input type="checkbox"/> Maxi+ 550 Low Seat height 22", handle 28.5-35"	<input type="checkbox"/> Maxi+ 650 Standard Seat height 25.5", handle 32-38"

Arm troughs or platform, bracket & pad adds 3" height to handle; does not include brake handles, choose brake below

Forearm trough (pad) One pad left One pad right 2 Pads, left & right side
If using pad without a brake, we will add a grip at end of pad

Platform, 27"w x 15.5"h Full platform, brakes mount at front edge, choose brake below

Brake handle position and style	One hand, operates brakes on both wheels	Braking both hands
Choose brake position	<input type="checkbox"/> Left mount one hand brake	<input type="checkbox"/> Right mount, one hand brake <input type="checkbox"/> Brakes on both sides
Choose style of brake action		<input type="checkbox"/> Joystick style, push/pull <input type="checkbox"/> Looped handle, squeeze/push

Alternate one hand braking systems, will not work with arm pads above

Crossbar, front push bar w/ one hand brake mounted center, replaces back strap, *not available for Maxi+* Left brake mount Right brake mount

HP Handle, one hand brake w/ angle adjustable bar on other side, *only available on Symphony, Left or right in prod # is based on HP side.* Left brake, right HP Right brake, left HP

Accessories	<input type="checkbox"/> D12533 Tray	<input type="checkbox"/> D12682 O2 Tank Holder	Ordering company _____ Return quote to _____ Email _____ Ph: _____ Zipcode for delivery _____
	<input type="checkbox"/> D13041 Slow Down Brakes	<input type="checkbox"/> P12519 Seat Pad,	
	<input type="checkbox"/> D12606 Cane holder	<input type="checkbox"/> D12523 Under seat bag	
	<input type="checkbox"/> D13235 Tall Handles for Legacy, adds 3" height		

Mark For _____ User weight _____ User height _____ Floor to Elbow _____

This walker will be built in our shop. It cannot be shipped totally assembled. The handles will need to be inserted and adjusted for height and length. It is not recommended for drop ship delivery unless a technician can be present to make the adjustments to fit the client. Arm troughs may block use of seat for sitting, platform will block use of seat. Note, the part numbers in the chart above do not indicate all components necessary to complete the walker assembly. Quote will show all necessary pieces. Non-returnable.

info@clarkehealthcare.com