



Custom Dolomite Rollator Worksheet for quote, return to info@clarkehealthcare.com or fax: 724-695-2922

1. Walker	Choose one model, size or color	4'10"-5'2" suggested height	5'2"-5'5" suggested height	5'4"-6'2" suggested height
	Legacy , <330 lbs cap. <div style="text-align: right;">Blue</div>	<input type="checkbox"/> Legacy 450 SuperLow Seat height 19", handle 28-31"	<input type="checkbox"/> Legacy 520 Low Seat height 22", handle 31-34"	<input type="checkbox"/> Legacy 600 Standard Seat height 24", handle 33-37"
	Symphony, <275 lbs cap. Silver frame	<input type="checkbox"/> Sym. 450 SuperLow Seat height 19", handle 25-30"	<input type="checkbox"/> Sym. 520 Low Seat height 22", handle 27.5-33"	<input type="checkbox"/> Sym. 600 Standard Seat height 24", handle 30-36"
	Maxi+, <452 lbs cap., Blue frame		<input type="checkbox"/> Maxi+ 550 Low Seat height 22", handle 28.5-35"	<input type="checkbox"/> Maxi+ 650 Standard Seat height 25.5", handle 32-38"

2. Arm troughs or platform , does not include brake handles, choose brake below	<input type="checkbox"/> One pad left	<input type="checkbox"/> One pad right	<input type="checkbox"/> 2 Pads, left & right side
Option 1 Forearm pad, adds 3" height	<i>If using pad without a brake, we will add a grip at end of pad</i>		
Option 2 Platform, adds 4" height	<input type="checkbox"/> Full platform, brakes mount at front edge, 27"w x 15.5"h		

3. Brake handle position and style	One hand, operates brakes on both wheels	Braking both hands
A. Choose brake position	<input type="checkbox"/> Left mount one hand brake	<input type="checkbox"/> Right mount, one hand brake
B. Choose style of brake handle and action	<input type="checkbox"/> Joystick style, push/pull	<input type="checkbox"/> Looped handle, squeeze/push
Alternate one hand braking systems, will not work with arm pads above		
Crossbar, front push bar w/ one hand brake mounted center, replaces back strap, <i>not available for Maxi+</i>	<input type="checkbox"/> Left brake mount	<input type="checkbox"/> Right brake mount
HP Handle, one hand brake w/ angle adjustable bar on other side, <i>only available on Symphony, Left or right in prod # is based on HP side.</i>	<input type="checkbox"/> Left brake, right HP	<input type="checkbox"/> Right brake, left HP

Accessories	D12533 Tray	D12682 O2 Tank Holder	Ordering company _____
	D13041 Slow Down Brakes	P12519 Seat Pad,	Return quote to _____
	D12606 Cane holder	D12523 Under seat bag	Email _____@_____
	D13235 Tall Handles for Legacy, adds 3" height		Ph: _____ Zipcode for delivery _____
	D13235M Tall Handles Maxi, adds 4" height		

Mark For _____ User weight _____ User height _____ Floor to Elbow _____

This walker will be built in our shop. It cannot be shipped totally assembled. The handles will need to be inserted and adjusted for height and length. It is not recommended for drop ship delivery unless a technician can be present to make the adjustments to fit the client. Arm troughs may block use of seat for sitting, platform will block use of seat. Note, the part numbers in the chart above do not indicate all components necessary to complete the walker assembly. Quote will show all necessary pieces. Non-returnable.