

Look for
Clarke Health Care Products
at upcoming CPTA
and APTA shows

Dear Physiotherapist

Now you can have it all -
Total Stairs & Gait Training solutions
in one perfect effective, space saving
machine!

You and your patients can now benefit from our latest development:
A combination of our well known Dynamic Stairs Trainer-DST,
together with Long Hand Rails-LHR and the incredible new
Convertible Stairs Ramp Curb-CSRC, gives you the ultimate tool for
stairs and gait training.

The new mini flip elevates the front ramp to provide an extended parallel bar walking distance. This configuration replaces regular parallel bars and saves valuable space in the ever-crowded physiotherapy room. Now you have an adjustable stairs trainer and parallel bars in one efficient machine for comprehensive treatment.



The height of the stairs is adjusted electronically to suit your patient's pace

The new mini flip ramp presents a practice curb while improving the floor-to-ramp transition. You can now choose between 3 options: moderate slope, three additional small stairs or long flat horizontal walking distance as parallel bars with curb practice.



Your distributor: 

Clarke Health Care Products, Inc.
7830 Steubenville Pike
Oakdale, PA 15071 USA
Ph: (888)347-4537 Fax: (724)695-2922
www.clarkehealthcare.com

Manufactured by:



DPE MEDICAL Ltd.
49 Yefe-nof, Shoeva 9085500 Israel
www.dpemed.com
info@dpemed.com

Scan the QR Codes with your smartphone
for amazing videos:

**Classic DST
in Action**



**CSRC
in Action**



**DST
for Children**



**The New DST
Compact**

