

Portable Grab Bars Retail 2022

Description	Overall length	Capacity	Part #	Retail
Shower head holder, short, adjustable angle	9"		R1400202S	\$131.00
Shower screen arm only, does not include curtain	28"-40.9"	N/A	R1400205S	\$510.00
Curtain only for R1400205S	46.8" x 31.2"	N/A	R14061208	\$146.00
Shower screen arm and curtain	11.7"-27.5"	N/A	R1499205S	\$543.00
One hand grip, fixed length	12.8"	231.5 lb	R1400220S	\$155.00
Small bar, telescoping	17"-21.4"	187.4 lb	R1400221S	\$249.00
Medium bar, telescoping	21.8"-26.7"	187.4 lb	R1400222S	\$261.00
Large bar, telescoping	25.7"-30.6"	187.4 lb	R1400223S	\$275.00
X-large bar, telescoping	34.5"-39.3"	176.4 lb	R1400224S	\$406.00
Small hand grip, fixed length	16.7"	220.5 lb	R1400225S	\$239.00
Medium hand grip, fixed length	22.2"	209.4 lb	R1400226S	\$239.00
Large hand grip, fixed length	33.9"	209.4 lb	R1400227S	\$239.00
QuatroPlus, 2 swivel points, 4 pads	14"-18.9"	352.8 lb	R1400241SU	\$444.00
QuatroPlus, 2 swivel points, 4 pads	18.7"-24.37"	352.8 lb	R1400242SU	\$485.00
QuatroPlus, 2 swivel points, 4 pads	22.6"-27.49"	352.8 lb	R1400243SU	\$563.00
QuatroPlus, 2 swivel points, 4 pads	21.4"-26.3"	352.8 lb	R1400244SU	\$613.00
QuatroPower tub, long, 3 swivel points, 4 pads	28.8"-38.2"	330.7 lb	R1400261SU	\$663.00
QuatroPower tub, long w/grip ext., 3 swiv pts, 4 pds	28.8"-38.2"	176 lb	R1400262SU	\$683.00
QuattroPower tub, short, 3 swivel pts, 4 pds	28"	176 lb	R1400263S	\$672.00
QuattroPower tub, short, w/grip ext. 3 swivel, 4 pds	33.15"	176 lb	R1400264S	\$677.00
QuatroPower support, 3 swivel pts, 4 pads	30.40"	176 lb	R1400265SU	\$677.00
QuatroPower support w/grip horn, 3 swiv pts, 4 pds	30.40"	176 lb	R1400266SU	\$697.00
Cardan joint adapter, pair, aluminum, swing and rotate possible, reduces capacity of bar to 110 lb.	3 5/8"	110 lb	R1400555	\$188.00
Diagonal joint adapter, pair, for Quatro Plus bars		264.6 lb	R1400557	\$188.00
Angle joint adapter, pair, hi-strength plastic, swings on one plane, reduces capacity to 143 lb.		165.3 lb	R1400558	\$119.00
Dual hook, white	3.5" drop	N/A	R1400581	\$26.00
Stainless steel plates, set 2, adhesive backing,	5.65"x5.65"		R1461414	\$52.00
Shower head holder swivel with indicator	11.7"-27.1"		R1499202S	\$206.00