

malte.

Smooth-running.

GMFCS
Levels III-IV

kg
max. 220



17.9"-48.3"

11.3"-44"



I am your versatile posterior walking

aid. I am always there for you - even from your first year of life. Walking is my biggest hobby. It doesn't matter where you want to go; I will follow you.

If you support yourself on my handles, you will quickly notice how I take on your body weight. In addition, my special abducted frame shape supports you so that you don't swing outwards when walking. My pelvic pelotte pad ensures that you walk upright and straight. Cool, aren't I?!





Buttocks pelotte pad
can become a backrest in
one single easy step.

Aluminium frame
can easily be folded
up.

Backwards stop
prevents inadvertent
rolling back.

EVA guide wheels
make the walking aid manoeuvrable.

Forearm rests
with handles.

Padded bench
can easily be folded
back.

Direction locks
facilitate
walking straight.



Abducted frame shape
supports you when walking
and prevents you from
swinging outwards.



Aluminium frame
can easily be folded
up via folding joints.

Frame colours
Size 0 mini



Frame colours
Size 1



Frame colours
Size 2



Frame colours
Size 3 and Size 4





Data	Size 0 mini*	Size 1	Size 2	Size 3	Size 4
Handle height	11.3"-20.7"	13.2"-22.6"	19.9"-29.2"	26.5"-36.6"	33.9"-44"
Height of forearm rests	-	17.9"-27.3"	24.5"-33.9"	30.8"-40.1"	38.2"-48.3"
Handle spacing	8.5"-14"	8.5"-14"	11.7"-17"	14.8"-20.2"	16.7"-22.2"
Seat height	-	12.8"	17.5"	23.4"	27.3"
Total width	20.2"	19.9"	23.8"	26.5"	29.2"
Total length	22.6"	22.6"	27.3"	30"	33.9"
Frame height (folded)	14.4"	9.7"	12.5"	13.2"	14.4"
Transport dimensions (L x W x H)	22.6"x20.2"x14.4"	22.6"x19.9"x9.75"	27.3"-23.8"x12.5"	30"x26.5"x13.2"	33.9"x29.2"x14.4"
max. load	44 lb.	66 lb.	110 lb.	165 lb./154 lb.	220 lb.
Weight	8.5lb.	15.6 lb.	18.3"	18.7 lb./21.6 lb.	22.4"
HMV No.	10.46.02.3029				

* Frame cannot be folded up

You can find more information on accessories in the Price List.



Distributed by Clarke Health Care Products 888-347-4537


schuchmann