

malte.

Smooth-running.

GMFCS  
Levels III-IV



### I am your versatile posterior walking

**aid.** I am always there for you - even from your first year of life. Walking is my biggest hobby. It doesn't matter where you want to go; I will follow you.

If you support yourself on my handles, you will quickly notice how I take on your body weight. In addition, my special abducted frame shape supports you so that you don't swing outwards when walking. My pelvic pelotte pad ensures that you walk upright and straight. Cool, aren't I?!





**Buttocks pelotte pad**  
can become a backrest in  
one single easy step.

**Aluminium frame**  
can easily be folded  
up.

**Backwards stop**  
prevents inadvertent  
rolling back.

**EVA guide wheels**  
make the walking aid manoeuvrable.

**Forearm rests**  
with handles.

**Padded bench**  
can easily be folded  
back.

**Direction locks**  
facilitate  
walking straight.



**Abducted frame shape**  
supports you when walking  
and prevents you from  
swinging outwards.



**Aluminium frame**  
can easily be folded up  
via folding joints.

**Frame colours**  
Size 0 mini



**Frame colours**  
Size 1



**Frame colours**  
Size 2



**Frame colours**  
Size 3 and Size 4







Data	Size 0 mini*	Size 1	Size 2	Size 3	Size 4
Handle height	11.3"-20.7"	13.2"-22.6"	19.9"-29.2"	26.5"-36.6"	33.9"-44"
Height of forearm rests	-	17.9"-27.3"	24.5"-33.9"	30.8"-40.1"	38.2"-48.3"
Handle spacing	8.5"-14"	8.5"-14"	11.7"-17"	14.8"-20.2"	16.7"-22.2"
Seat height	-	12.8"	17.5"	23.4"	27.3"
Total width	20.2"	19.9"	23.8"	26.5"	29.2"
Total length	22.6"	22.6"	27.3"	30"	33.9"
Frame height (folded)	14.4"	9.7"	12.5"	13.2"	14.4"
Transport dimensions (L x W x H)	22.6"x20.2"x14.4"	22.6"x19.9"x9.75"	27.3"-23.8"x12.5"	30"x26.5"X13.2:	33.9"x29.2"x14.4"
max. load	44 lb.	66 lb.	110 lb.	165 lb./154 lb.	220 lb.
Weight	8.5lb.	15.6 lb.	18.3"	18.7 lb./21.6 lb.	22.4"
HMV No.	10.46.02.3029				

\* Frame cannot be folded up



You can find more information on accessories in the Price List.

Distributed by Clarke Health Care Products

888-347-4537

