

marcy.

There for you from the start.

GMFCS
Levels II-IV

kg
max. 65



41-98 cm

30-88 cm

Handles
can be adapted in height
and width.

Storage basket
provides space for toys or
drinks.

Aluminium frame
highly compact; sizes 2+3 can be folded.


PU tires
are highly resistant.

Friction brakes on the rear wheels throttle the speed so that
you can concentrate completely on walking.

Direction locks
support you
in walking straight ahead.

I am your anterior walking aid, and will assist you in your very first attempts to walk.

I can hardly wait to start off with you! I will support your balance and simultaneously train your torso musculature. You can adapt my wheelbase to suit: if you adjust it so that it is somewhat shorter, I become highly manoeuvrable and you can move me better. I give you more security and support with a long wheelbase. If you find it hard to walk in a straight line, I am happy to assist you with my direction locks.





Three dimensionally-adjustable forearm supports
stabilise your gait pattern.



Optional parking handbrake
provide you with more safety and control – also in combination with the forearm supports.



Data	Size 0	Size 1	Size 2	Size 3
Handle height	11.7"-18.3"	14.8"-23"	20.6"-29.25"	24.5"-34.3"
Height of forearm rests	15.9"-22.6"	19.1"-27.3"	24.5"-33.1"	28.4"-38.2"
Handle spacing	7"-12.4"	8.5"-14"	11.7"-17.1"	14"-19.5"
Seat height	9.75"	11.7"	17.5"	17.5"
Total width	17.1"	18.7"	24.1"	27.6"
Total length	14.2"	16.5"	22.4"	23.7"
Total length (telescopic wheelbase)	15.2"-18.3"	17.5"-21.4"	23.4"-28.8"	25.3"-30.8"
max. load	44 lb.	66.1 lb.	110 lb.	143.3 lb.
Weight	9.2 lb.	11 lb.	16.3 lb.	20 lb.
HMV No.	10.46.02.3049			

Frame colours
Size 0

- Telemagenta
- Sky blue

Frame colours
Size 1

- Telemagenta
- Sky blue

Frame colours
Size 2

- Lazurite blue metallic
- Traffic red

Frame colours
Size 3

- Anthracite metallic
- Pure white

You can find more information on accessories in the Price List.



Clarke Health Care Products 888-347-4537

