



GMFCS  
Levels II-IV



**I am your anterior walking aid, and will assist you in your very first attempts to walk.**

I can hardly wait to start off with you! I will support your balance and simultaneously train your torso musculature. You can adapt my wheelbase to suit: if you adjust it so that it is somewhat shorter, I become highly manoeuvrable and you can move me better. I give you more security and support with a long wheelbase. If you find it hard to walk in a straight line, I am happy to assist you with my direction locks.





**Handles**

can be adapted in height and width.

**Aluminium frame**

highly compact; sizes 2+3 can be folded.

**Friction brakes on the rear wheels**

throttle the speed so that you can concentrate completely on walking.

**Storage basket**

provides space for toys or drinks.

**Direction locks**

support you walking straight ahead.

**PU tires**

are highly resistant.



**Three dimensionally-adjustable forearm supports**  
stabilise your gait pattern.

**Optional parking handbrake**

provide you with more safety and control - also in combination with the forearm supports.

**Frame colours  
Size 0****Frame colours  
Size 1****Frame colours  
Size 2****Frame colours  
Size 3**



Data	Size 0	Size 1	Size 2	Size 3
Handle height	11.7"-18.3"	14.8"-23"	20.6"-29.25"	24.5"-34.3"
Height of forearm rests	15.9"-22.6"	19.1"-27.3"	24.5"-33.1"	28.4"-38.2"
Handle spacing	7"-12.4"	8.5"-14"	11.7"-17.1"	14"-19.5"
Seat height	9.75"	11.7"	17.5"	17.5"
Total width	17.1"	18.7"	24.1"	27.6"
Total length	14.2"	16.5"	22.4"	23.7"
Total length (telescopic wheelbase)	15.2"-18.3"	17.5"-21.4"	23.4"-28.8"	25.3"-30.8"
max. load	44 lb.	66.1 lb.	110 lb.	143.3 lb.
Weight	9.2 lb.	11 lb.	16.3 lb.	20 lb.
HMV No.	10.46.02.3049			



You can find more information on accessories in the Price List.



Clarke Health Care Products 888-347-4537

