



The MAXI+

Maxi+ is designed for the heavier, larger, or taller user. The wider frame provides more seat room and a larger walking area. Dolomite has the strongest frames on the market; Maxi+ is rated for a 452 lb. capacity. Choose Maxi+ for a stable walking experience and travel safely in style.

Components

Maxi+ is made from high tensile steel with welded joints for strength and low maintenance. Maxi+ comes in Metallic Blue epoxy-based powder coat finish. Dolomite Maxi+ comes in two frame heights; see chart for specs. The frames fold front to back for storage. A safety latch prevents accidental folding when using the seat handle.

Maxi+ includes a wide flip-up seat for more walking room. Anatomical grips have a larger surface, which reduces hand fatigue for arthritic users. Flexible back strap is adjustable. Large 8" wheels roll easily indoors or outdoors.

Braking

Maxi+ comes with the patented DABS braking system, with adjustable tension for varied grip strength. Brake handles have a wide flange and large loops for comfortable hand positioning. The brake mechanism is located within the fork for protection. When parked and using the seat, push down the handles to set the parking brakes, audible click is heard.

Features

- Extra wide, 20.5" between handles
- Large flip-up seat
- 452 pound capacity
- Perfect for the taller person, 25" seat height
- Large 8" wheels
- Handles adjust 7" in 1" increments
- Adjustable back support

Warranty

To original purchaser, Dolomite provides a limited lifetime warranty against manufacturer's defects on the frame, 5 years on brake handle and one year on parts. Tires, brake cables and grips are excluded from this warranty.



CLARKE HEALTH CARE PRODUCTS, INC.

7830 Steubenville Pike • Oakdale, PA 15071

Phone (888) 347-4537 • Fax (724) 695-2922 • www.clarkehealthcare.com



The MAXI+

Accessories



TRAY



OXYGEN TANK HOLDER



CANE HOLDER



TALL HANDLES
MAXI+



SLOW DOWN BRAKE



ONE HAND BRAKE



UNDER SEAT BAG



SEAT AND BACK PAD

Accessories can enhance your walker – the Tray can carry food, Oxygen Tank Holder will safely support a B, C or D tank. Cane Holder can secure a care of crutch on the side of the walker. Slow Down Brakes have a tension control to allow individual settings for each wheel. One-Hand Brakes allow braking from one handle. Under Seat Bag gives you a place to stash items out-of-sight. Handle extensions add 4" height. Large basket for personal items.

Specifications:

| SPECIFICATIONS | LOW | STANDARD |
|---------------------------|---------------------------------------|-----------------|
| Dist. Between Rear Wheels | 22" | 22" |
| Handle Height | 28.5" - 35" | 32.5" - 39.5" |
| Tall Handle Height | 32.5" - 39" | 36" - 43" |
| Seat Height from Floor | 22" | 25" |
| Seat Dimensions | 8.5" Deep x 17" Wide | |
| Overall W x L | 27.5" W x 26" L | 27.5" W x 29" L |
| Width Between Handles | 20.5" | 20.5" |
| Folded Dimensions | 33" Length x 27.5" Width x 11" Height | |
| Walker Weight (lbs.) | 20.5 | 21 |
| Weight Capacity (lbs.) | 452 | 452 |
| Product Number | D12122 | D12120 |

CLARKE HEALTH CARE PRODUCTS, INC.

7830 Steubenville Pike • Oakdale, PA 15071

Phone (888) 347-4537 • Fax (724) 695-2922 • www.clarkehealthcare.com

