



# The MAXI+

Maxi+ is designed for the heavier, larger, or taller user. The wider frame provides more seat room and a larger walking area. Dolomite has the strongest frames on the market; Maxi+ is rated for a 452 lb. capacity. Choose Maxi+ for a stable walking experience and travel safely in style.

## Components

Maxi+ is made from high tensile steel with welded joints for strength and low maintenance. Maxi+ comes in Metallic Blue epoxy-based powder coat finish. Dolomite Maxi+ comes in two frame heights; see chart for specs. The frames fold front to back for storage. A safety latch prevents accidental folding when using the seat handle.

Maxi+ includes a wide flip-up seat for more walking room. Anatomical grips have a larger surface, which reduces hand fatigue for arthritic users. Flexible back strap is adjustable. Large 8" wheels roll easily indoors or outdoors.

## Braking

Maxi+ comes with the patented DABS braking system, with adjustable tension for varied grip strength. Brake handles have a wide flange and large loops for comfortable hand positioning. The brake mechanism is located within the fork for protection. When parked and using the seat, push down the handles to set the parking brakes.

## Features

- Extra wide, 20.5" between handles
- Large flip-up seat
- 452 pound capacity
- Perfect for the taller person, 24" seat height
- Large 8" wheels
- Handles adjust 7" in 1" increments, up to 39.5" without need of additional tubes
- Large basket for personal items
- Adjustable back support

## Warranty

To original purchaser, Dolomite provides a lifetime warranty against manufacturer's defects on the frame, 5 years on brake handle and one year on parts. Tires, brake cables and grips are excluded from this warranty.



**CLARKE HEALTH CARE PRODUCTS, INC.**

7830 Steubenville Pike • Oakdale, PA 15071

Phone (888) 347-4537 • Fax (724) 695-2922 • [www.clarkehealthcare.com](http://www.clarkehealthcare.com)



# The MAXI+

## Accessories



TRAY



OXYGEN TANK HOLDER



CANE HOLDER



SLOW DOWN BRAKE



ONE HAND BRAKE



UNDER SEAT BAG

Accessories can enhance the use of your walker – the Tray gives you a platform to carry food, the Oxygen Tank Holder will safely support a B, C or D tank. A Cane Holder can secure a cane or crutch to the side of the walker, Slow Down Brakes have a tension control to allow individual settings for each wheel. One-Hand Brakes can control your braking from one handle. Under Seat Bag gives you a place to stash items out of sight.

### Specifications:

SPECIFICATIONS	LOW	STANDARD
Dist. Between Rear Wheels	22"	22"
Minimum Handle Height	28.5"	32.5"
Maximum Handle Height	35"	39.5"
Seat Height from Floor	22"	24"
Seat Dimensions	8.5" Deep x 17" Wide	
Overall Width	27.5"	27.5"
Overall Length	26"	29"
Width Between Handles	20.5"	20.5"
Folded Dimensions	33" Length x 27.5" Width x 11" Height	
Walker Weight (lbs.)	20.5	21
Weight Capacity (lbs.)	452	452
Product Number	D12122	D12120

### CLARKE HEALTH CARE PRODUCTS, INC.

7830 Steubenville Pike • Oakdale, PA 15071

Phone (888) 347-4537 • Fax (724) 695-2922 • [www.clarkehealthcare.com](http://www.clarkehealthcare.com)

