

# AQUATEC® *The* OCEAN E-VIP



Who is taking care of the caregiver? Aquatec is, with the Ocean E-VIP providing seating and positioning in an adjustable shower commode chair system, designed to alleviate lifting and back strains. Ocean E-VIP consists of a seating system mounted on an electronically height-adjustable wheeled chassis.

For all your hygiene needs, Ocean E-VIP rolls over a commode or into a shower. The seat rises to 41" with a remote control. With the same ease, the remote controls the tilt angle of the seat. The patient is placed at a height or angle best suited for the care necessary. The chair can operate a full day, up to 40 cycles on a single charge. There are plenty of positioning adjustments made with seat, back, head and foot rests. All this movement is done with minimal effort from the caregiver. Ocean E-VIP is suitable for home or institutional use, anywhere a caregiver needs more assistance.

## Features

- Adjust seat height by remote
- Tilt-in-space position controlled by remote
- Corrosion-free materials
- Adjustable head support and footrests
- Padded chest belt included
- Directional lock on rear wheel for better handling
- Complete chair can be cleaned in autoclave
- Rechargeable battery provides 40 cycles on a single charge

## Components

The frame is constructed from stainless steel, aluminum or plastics to achieve a corrosion free product. The electrical components are sealed, allowing the chair to be cleaned or disinfected by commercial cleaning methods. Recharging can be done on the chair, or when used in shower rooms, the battery case can be removed and recharged off the chair.

Three wheels have foot activated brakes and one rear wheel has a directional lock for better handling when transporting the patient. Padded chest belt is included. The head support can be adjusted horizontally and vertically. Removable, washable back upholstery is also adjustable for best fit.



*Sealed electrical components*

**CLARKE HEALTH CARE PRODUCTS, INC.**

1003 International Drive • Oakdale, PA 15071

Phone (888) 347-4537 • Fax (724) 695-2922 • [www.clarkehealthcare.com](http://www.clarkehealthcare.com)

# *The* OCEAN E-VIP Features



*Lowest Height*



*Raised, Tilted*



*Seat Raised*



*Full Height, Tilted*



*Full Height*



*Padded Chest Belt*



*Adjustable Footplate*



*Directional Caster Lock*



*Adjustable Head Support*

**AQUATEC®**



# The OCEAN E-VIP Accessories



Neckrest Pad



Headrest Pad



Side Lateral with Height Adjustable Pad



Swing Away Side Lateral



Side Lateral with Angle Adjustable Pad



Arm Trough 16"



Arm Trough 11"



Arm Width Extension



Special Soft Seat Overlay



Toilet Overlay



Action Gel Overlay



Splash Guard



Solid Insert



Padded Seat Belt



Safety Harness Combination



Collection Pan with Lid

## Warranty

Ocean E-VIP has a limited two year warranty against manufacturer's defects. The battery is excluded.

SPECIFICATIONS	Ocean E-VIP
Overall height, without headrest	41.5" - 63"
Overall width	28.5"
Overall length including footrest	39.5"
Seat height	19.5" - 41"
Seat width	18.75"
Seat depth	17.5"
Width between casters	17.5"
Width between armrests	18"
Width between armrests with opt. arm width extension	21"
Battery	24V, 2.9Ah
Weight of chair	83.75 lb.
Capacity	330 lb.
Product number	A1471208



Safety Support Bar



XL Safety Bar



Arm Lock



Calf Support



Amputee Support



Calf Strap

## CLARKE HEALTH CARE PRODUCTS, INC.

1003 International Drive • Oakdale, PA 15071

Phone (888) 347-4537 • Fax (724) 695-2922 • [www.clarkehealthcare.com](http://www.clarkehealthcare.com)

