

**AQUATEC®**

# The OCEAN VIP / OCEAN DUAL



## *The OceanVIP or OceanDual with Recline Back*

The next wave in comfort and style, Aquatec's OceanVIP tilting shower chair is height adjustable, for easier transfers or toilet clearance. OceanVIP provides correct posture for your hygiene needs with tilt-in-space or recline positioning choices. Recline positioning now available on OceanDual or with the optional back support for your existing OceanVIP.

### *Features*

- Stainless steel frame
- Seat height easily adjusted from 19.5" to 25.5" from floor
- Adjustable headrest
- Adjustable, machine washable back upholstery
- 5" casters with locks
- Height adjustable, flip-up armrests

### *Components*

OceanVIP is built on a stainless-steel frame with corrosion-free components. The seat is waterfall style with a hygiene cut-out for comfort. Soft overlays give you seating choices for toileting or showering. Two stainless-steel gas springs will lift 330 pounds. Adjustable seat height is set with easy-to-move clips. Footrests can be lengthened easily with locking pins and include tilting footplate and heel loops. The back upholstery has 5 buckles to adjust tension for seating comfort.

OceanDual, with reclining back support, provides even greater seating and positioning possibilities. The angle between the reclining backrest and seat is infinitely adjustable from 10° to 50° for a total angle of 85° for the complete chair. The operating lever locks to avoid unintentional angle movement. A pair of anti-tippers is included with this back option.

The OceanVIP rolls smoothly on 5" casters. Three wheels have a foot activated brake and one rear wheel has direction lock for straight line travel. A handle on the back push bar sets the tilt from 0°-35°. Armrests flip-up for side transfers and can be positioned at two heights. Optional arm locks are available.

### *Warranty*

OceanVIP shower commode chair has a limited two year warranty against manufacturer's defects, excluding soft seat and upholstery. The optional reclining back support has a limited two year warranty against manufacturer's defects. Cylinders are warranted for one year.



*OceanDual in full tilt and recline, showing optional amp support*



*OceanVIP in full tilt*



*OceanDual*

## The OCEAN VIP / OCEAN DUAL Features



*Tilt Angle Adjustment*



*Seat Height Adjustment*



*Adjustable Upholstery*



*Adjustable Armrest Height*



*Head Support*



*Collection Pan with Lid*



*Directional Lock*



*Adjustable Footrest and Heel Loops*

**AQUATEC®**

**CLARKE HEALTH CARE PRODUCTS, INC.**

1003 International Drive • Oakdale, PA 15071

Phone (888) 347-4537 • Fax (724) 695-2922 • [www.clarkehealthcare.com](http://www.clarkehealthcare.com)

# The OCEAN VIP / OCEAN DUAL Accessories



Toilet Overlay



Action Gel Overlay



Special Soft Seat Overlay



Solid Insert



Splash Guard



Calf Strap



Safety Harness Combination



Padded Seat and Chest Belt



Safety Support Bar



XL Safety Bar



Calf Support



Amputee Support



Swing Away Side Lateral



Side Lateral with Angle Adjustable Pad



Side Lateral with Height Adjustable Pad



Arm Width Extension



Headrest Pad



Neckrest Pad



Arm Lock



Arm Trough 16"



Arm Trough 11"

## SPECIFICATIONS

	OceanVIP
Top of back to seat	21"
Seat height from floor	19.5" / 21.5" / 25.75"
Seat area (L xW)	17.5" x 19.5"
Overall width of chair	22.5"
Overall length 0° tilt	41.25"
Overall length, full tilt, footrest at max, headrest at max	60.2"
Width between arms	17.5"
Height floor to lower seat edge w/o pan holder	17.25" - 22.25"
Tilt angle	0°-35°
Distance between rear wheels	18"
Max weight capacity	330 lbs.
Weight of chair	48 lbs.
Product Number (Tilt Chair)	A1470713

## RECLINING BACK SPECIFICATIONS

Width between arms at end of handgrip	18.5"
Reclining back angle	0°-50°
Overall width of back	22.5"
Overall length: chair with recline back, no tilt, no recline, from anti-tipper to footrest	42"
Overall length: full tilt, no recline push bar to foot plate (no head support)	54"
Overall length: full tilt and recline, head support to foot plate	68"
Height from seat to top of upholstery	20"
Product Number (Recline Back only)	A1470085
Product Number (OceanDual - Tilt and Recline)	A1470829

## CLARKE HEALTH CARE PRODUCTS, INC.

1003 International Drive • Oakdale, PA 15071  
 Phone (888) 347-4537 • Fax (724) 695-2922 • www.clarkehealthcare.com

The OCEAN VIP / OCEAN DUAL