

AQUATEC® *The* OCEAN

Aquatec Ocean has been given a facelift, enhanced with features to better fit your needs.

Features

- Push bar has been moved upward so the chair can get closer to the wall when positioned over the toilet.
- Backrest frame has been redesigned so that the same soft seat overlays can be used across the entire Ocean line.
- Backrest upholstery straps have been lengthened to allow greater adjustment for comfort and positioning.
- Increased footplate dimension; is now angle adjustable.
- A common armrest with bracket is now in place on all Oceans. The same arm locks can now be used on Ocean and OceanVIP.
- The XL model now has a wider backrest frame for more comfort.
- Height adjustable seat; no tools needed.

Components

Ocean is a family of stainless steel, height adjustable shower/commode chairs. Changing the seat height is done with push clips in 1" increments for a total of 5" height range, no tools required. Seat frame has a 10° static tilt. All models have 5" casters with wheel locks. PVC-coated, nylon backrest upholstery is washable and tension is easily adjusted. Coated urethane seat overlay has a keyhole shape cutout. Armrests flip up for transfers and can be adjusted upwards by 1 1/2". Swing-away footrests include heel loops. They can be adjusted in length from seat and footplate angle. A collection pan with lid is included. Entire chair can be disassembled for cleaning in a tub or autoclave.

OceanSP is a self-propelled model with 24" wheels, hand brakes and push rims. The large wheels have quick-disconnect pins, allowing you to remove the wheels and use the 5" casters. OceanXL has a wider backrest, total of 21" between arms, for greater comfort. The arm brackets are reinforced for greater weight capacity.

Warranty

Aquatec Ocean shower/commode chair frame is warranted for two years against manufacturer's defects, excluding upholstery and soft seat overlay.

CLARKE HEALTH CARE PRODUCTS, INC.

1003 International Drive • Oakdale, PA 15071
 Phone (888) 347-4537 • Fax (724) 695-2922 • www.clarkehealthcare.com



OceanSP



OceanXL

The OCEAN Features



Adjustable Footrest and Heel Loops



Adjustable Back Upholstery



Adjustable Seat Height



Adjustable Armrest Height



Seat Overlay with Hygiene Cutout



Collection Pan with Lid



Disassembles for Travel

AQUATEC®

The OCEAN Accessories



Toilet Overlay



Action Gel Overlay



Special Soft Seat Overlay



Key Insert



Swing Away Side Lateral



Side Lateral with Angle Adjustable Pad



Side Lateral with Height Adjustable Pad



Arm Lock



Splash Guard



Calf Strap



Safety Harness Combination



Padded Seat and Chest Belt



Carry Bag



Safety Support Bar



Safety Bar for XL



Amputee Support



Calf Support



Arm Trough 16"



Arm Trough 11"

SPECIFICATIONS	Ocean (5" casters)	OceanSP (24" wheels)	OceanXL
Width of chair	22.5"	27"	24.5"
Width between arms	18"	18"	21"
Height of arm, from seat	8.5" / 10"	8.5" / 10"	8.5" / 10"
Length of chair	34"	42.5"	34"
Height of chair	38" - 42.5"	38" - 42.5"	38" - 42.5"
Seat area (L x W)	19" x 18"	19" x 18"	19" x 18"
Seat height from floor	19" - 23"	19" - 23"	19" x 23"
Standard wheel diameter	5"	5" + 24"	5"
Seat height adjustment increments	1"	1"	1"
Height under seat	16" - 21"	16" - 21"	16" - 21"
Min. width between side sections, under frame, for toilet clearance	15"	15"	15"
Capacity	290 lb.	290 lb.	400 lb.
Weight of chair	33.5 lb.	47 lb.	33.5 lb.
Back support (L x W)	14" x 16"	14" x 16"	14" x 21"
Standard hygiene seat opening (L x W)	15.5" x 8"	15.5" x 8"	15.5" x 8"

CLARKE HEALTH CARE PRODUCTS, INC.

1003 International Drive • Oakdale, PA 15071
 Phone (888) 347-4537 • Fax (724) 695-2922 • www.clarkehealthcare.com

The OCEAN