

Aquatec® | *Shower commode chairs*



Aquatec ***Ocean, Ocean XL, Ocean Vip, Ocean Dual Vip, Ocean E-Vip***

The Aquatec ***Ocean*** series of shower commode chairs offers the right model for all demands.

AQUATEC®

 **Clarke**
Health Care Products

www.clarkehealthcare.com
888-347-4537



Aquatec **Ocean**

- Stable, stainless steel frame
- Individually adjustable seat height: no tools required
- Machine-washable adjustable backrest with quick-release buckles, can be adjusted to the user
- Space-saving push bar to enable optimal positioning over the toilet
- Ergonomic, swing-away armrests can be adjusted to two different height levels
- Height-adjustable, swing-away and detachable footrests for simple transfer
- Fold-away, textured footplates with heel loops provide a safe position

The **Ocean** is available as follows:

- Standard version with castors
- Self-propel model with 24" wheels, quick-release-axles and four castors
- XL version (400 lb. load capacity and more space between the armrests) with castors
- All models can be retrofitted with 24" wheels and/or to the XL version

Options and features



Ocean 24" wheels
Self-propel model with 24" wheels, quick-release axles and four castors.



Ocean XL
More space between the armrests with a load capacity of up to 400 lb..



Individually adjustable seat height: no tools required.



Easy disassembly for mobile use, e.g. on travels.



Stable, stainless steel frame.



Aquatec **Ocean Vip**

- Stable, stainless steel frame
- Seat tilt can be individually adjusted (tilt angle 0° - 35°)
- Stainless steel gas spring
- Extremely stable even at maximum tilt angle
- Seat height adjustable in 3 steps, no tools required
- Very comfortable, ergonomically contoured seat
- Machine-washable adjustable backrest with quick-release buckles – can be easily adjusted to the user

The **Ocean Vip** is available as follows:

- Standard version with castors
- Standard version with castors and toilet pan holder
- Standard version with castors, toilet pan holder and soft seat
- The **Ocean Vip** can be retrofitted with XL armrests.

Options and features



Seat tilt

Tilt angle 0° - 35°.



Headrest

Height and depth-adjustable headrests included.



Directional lock

Easy to operate brakes on 3 castors, 1 castor with directional lock to help push the chairs down straight floors.



Adjustable seat height

Adjustable in 3 steps, no tools required.



Aquatec **Ocean Dual-Vip**

- Seat tilt can be individually adjusted (tilt angle 0° - 35°)
- Two additional stainless steel gas springs to lower the backrest (recline angle 10° - 50°)
- Second operating lever with safety lock to lower the backrest
- Particularly stable even at maximum tilt angle, anti-tippers included
- Machine-washable adjustable backrest with quick-release buckles – can be adjusted to the user

The **Ocean Dual Vip** is available as follows:

- Standard version with castors
- The **Ocean Vip** can be retrofitted with reclining backrests for the **Ocean Dual Vip** and be fitted with anti-tippers
- The **Ocean Dual Vip** can be retrofitted with XL armrests

Options and features



Seat surface 0°, backrest 25°



Full seat tilt 35°, back rest 50°



Anti-tippers included



Aquatec **Ocean E-Vip**

- Seat height can be steplessly adjusted between 500 and 1050 mm electrically
- Seat tilt can be electrically tilted forwards and backwards (tilt angle -5° to 35°)
- Particularly stable even at maximum tilt angle and maximum seat height
- Clearly laid out hand control for simple adjustment of the seat height and seat tilt (4 buttons only)
- Sealed electronic components: ideal for use in the shower
- Can be cleaned in automatic cleaning stations at temperatures up to 85 °C!
- Detachable battery pack can be recharged easily

The **Ocean E-Vip** is available as follows:

- Standard version with castors
- The **Ocean E-Vip** can be retrofitted with XL armrests

Options and features



Padded chest belt included



Maximum tilt angle and maximum seat height



Detachable battery pack



Clearly laid out hand control

Accessories Aquatec Ocean



Soft seat

With key-shaped cut out.
A16341 for **Ocean / Ocean XL / Ocean E-Vip**
A1470075 for **Ocean Vip / Ocean Dual Vip**



Soft seat insert

For soft seat with key-shaped cut out.
A16343 for all models of the **Ocean** series



Soft seat

With small, oval opening. **A16342** for all models of the **Ocean** series



Splash guard

Simply snap-on, for soft seat with key-shaped cut out.
A16511 for all models of the **Ocean** series

Accessories Aquatec Ocean



Universal soft seat
Closed seat surface.
A1526032 for all models of the **Ocean** series



Ergonomic special soft seat
Extra soft, textured surface
A1531259 soft for **Ocean**, **Ocean XL**, **Ocean Vip**, **Ocean Dual Vip**



Variable soft seat
Hygiene cut out can be placed in all 4 directions. Easy and tool less assembly.
A1544469 for **Ocean/Ocean Vip / Ocean Dual Vip**



Safety bar
For additional safety. The surface with soft padding. Easy to assemble.
D10167 for all models of the **Ocean** series. Length: 540 mm
A1470078 XL safety bar for all XL versions of the **Ocean** series. Length: 620 mm



Safety belts
For safe positioning. Individually adjustable, high load capacity and easy to assemble. Machine-washable. For all models of the **Ocean** series. **A1470081** Chest belt, grey **A1470082** Pelvic belt, grey



Armrest locking device
To lock the armrest.
A1470077 for all models of the **Ocean** series



XL back
Increases the distance between the armrests by 3". Retrofit kit for **Ocean**.
A1534979



Backrest
Solid backrest with ergonomic shape, easy to clean and dry. Snap-on assembly.
A1555323 for **Ocean/Ocean Vip / Ocean Dual Vip** and **Ocean E-Vip**



Soft backrest cushion
For additional stability and comfort - also for small users and young adults
A13681 for **Ocean / Ocean XL / Ocean Vip / Ocean Dual Vip** and **Ocean E-Vip**



Wider XL armrests
Increases the distance between the armrests by 1-1/2" on each side.
Retrofit kit for **Ocean Vip**, **Ocean Dual Vip** and **Ocean E-Vip**
A1470086 Wider XL armrests (pair)



24" wheels (set)
Incl. brakes and removable axles.
A16379 for **Ocean**



Calf-and amputee support
Padded, both height and angle adjustable.
A1556894 left for **Ocean / Ocean XL / Ocean Vip / Ocean Dual Vip** and **Ocean E-Vip**
A1556893 right for **Ocean / Ocean XL / Ocean Vip / Ocean Dual Vip** and **Ocean E-Vip**

Accessories Aquatec Ocean



Calf support with foot support
Padded, both height and angle adjustable with additional foot support.

A1546705 left for *Ocean / Ocean XL/Ocean Vip/Ocean Dual Vip* and *Ocean E-Vip*

A1546439 right for *Ocean / Ocean XL/Ocean Vip / Ocean Dual Vip* and *Ocean E-Vip*



Lateral
For safe positioning. Width and height adjustable – also for small users and young adults.

A1535077 for *Ocean / Ocean XL / Ocean Vip/ Ocean Dual Vip* and *Ocean E-Vip*



Ergonomic head rest holder
Safe and comfort with various adjustment positions (pad not included).

A1531613 for *Ocean Vip / Ocean Dual Vip* and *Ocean E-Vip*



Head/Neckrest (neckrest shown)

A1531465 Headrest for *Ocean Vip, Ocean Dual Vip* and *Ocean E-Vip*

A1531608 Neckrest for *Ocean Vip, Ocean Dual Vip* and *Ocean E-Vip*



Front anti tipper

For additional safety

A1535936 for *Ocean Vip / Ocean Vip XL / Ocean Dual Vip*



Side Lateral

Pad 7" x 4.3", extends 5.8"
O434F3 for all models of Ocean series.



Side Lateral

Swing away

RP05261197 fits all models of Oceans



Arm Trough

Set of 2 pads with hardware
18" long

O436A3 fits all models of Oceans



Action Gel Overlay

Gel pad for additional seat cushioning. 18x18 with cutout

D100031 for all models of Oceans



Calf Strap

4" wide strap, secures with Velcro.

D10175 for all models of Oceans

Technical data

Order Numbers	Aquatec Ocean	Aquatec Ocean 24" Wheels	Aquatec Ocean XL	Aquatec Ocean Vip With Tilt	Aquatec Ocean Dual Vip With Tilt & Recline	Aquatec Ocean E-Vip
Shower Commode Chair With Castors						A1471208
Shower Commode Chair With 24" Wheels		A1534327SP				
Shower Commode Chair With Castors, Toilet Pan Holder and Soft Seat	A1534327		A1534328	A1470713	A1525707	

Dimensions						
Seat Width (A)	18.7"	18.7"	18.7"	19.5"	19.5"	18.7"
Seat Depth (B)	17.5"	17.5"	17.5"	17.5"	17.5"	17.5"
Seat Height (C)	-	-	-	19.9" - 23.8"	19.9" - 23.8"	19.5" - 41"
Total Width (D)	21.8"	26.3"	25"	23.2"	24.3"	28.5"
Total Depth (E)	35"	41.3"	35"	41.3"	41.3"	39.4"
Total Height without Adjustable Backrest (F)	36.6" - 41.7"	36.6" - 41.7"	36.6" - 41.7"	42.3" - 46.2"	42.3" - 46.2"	41.7" - 63"
Depth without Footrests (G)	27.3"	33.3"	27.3"	33"	33"	39.4"
Footplate Depth (H)	6"	6"	6"	6.2"	6.2"	6"
Height to Armrest (I)	26.5" - 33"	26.5" - 33"	26.5" - 33"	28" - 33.5"	28" - 33.5"	28.5" - 51.5"
Height, Footrest to Seat (J)	10" - 19.3"	10" - 19.3"	10" - 19.3"	13.5" - 18"	13.5" - 18"	13.5" - 18"
Height, Ground to Lower Seat Edge With Holder for Toilet Pan (K)	15.2" - 20.3"	15.2" - 20.3"	15.2" - 20.3"	16.4" - 20.3"	16.4" - 20.3"	16.4" - 37.8"
Without Holder for Toilet Pan (K)	16" - 21"	16" - 21"	16" - 21"	17" - 21"	17" - 21"	18" - 39.6"
Height, Seat to Armrest (L)	8.2" / 9.75"	8.2" / 9.75"	8.2" / 9.75"	8.4" / 9.75"	8.4" / 9.75"	9" - 10.5"
Width between Armrests (M)	18"	18"	21"	18.5"	18.5"	18.3"
Width Cutout in Seat Plate	8"	8"	8"	8"	8"	8"
Depth Cutout in Seat Plate	12.5"	12.5"	12.5"	12.5"	12.5"	12.5"
Tilt Angle Seat	-	-	-	0° - 35°	0° - 35°	-
Tilt Angle Back	-	-	-	-	10° - 50°	-
Graduation Height Adjustment	1"	1"	1"	2"	2"	stepless
Diameter Castors	5"	5"	5"	5"	5"	5"
Diameter Self Propel Drive Wheels	-	24"	-	-	-	-
Width Between 24" Wheels (N)	-	21"	-	-	-	-
Weight	32 lbs.	43 lbs.	32 lbs.	52 lbs.	62.4 lbs.	86 lbs.
Load Capacity	330 lbs.	330 lbs.	330 lbs.	330 lbs.	330 lbs.	330 lbs.
Color	White	White	White	White	White	White

Heel loops are included in the scope of delivery.

