

DOLOMITE | *STEP UP*



The rehab optimizer

DOLOMITE STEP UP is a new, modern platform support. Made to provide rehabilitation to several patients every day. Adjust to optimal ergonomic position with a few simple steps – quick and easy! **DOLOMITE STEP UP** can be manually adjusted or hydraulically. Assembly of accessories is tool free. Are you ready to step up?





Made for training

The base frame is designed to enhance premium gait area. The raised frame in front allows longer step. The width inside the frame is generous.

Ergonomic adjustment of arm support

Arm pads are in soft material for pressure relief. Securely assembled with microfilm. Quick and simple to remove for cleaning or replacement.

Transfer and access

Well thought out for transfer situations. From bed, the lowered rear leg helps backing under the bed.

When transferring from wheelchairs benefit from the generous frame size and the width adjustment.

Features and Options



By releasing the central knob, both arm pads are simultaneously adjusted in width. The position of the arm supports remains at an angle of 25 degrees from the shoulder points - the ideal position for keeping upper strength.

Accessories and options

Quick and easy assembly in track or prepared recess.



Practical basket
Keeps belongings within easy to reach range.
Part no D1531390



Tray
For small things .Fits on basket.
Part no D1531391



Elbow stop
Secures elbows from falling out of position.
Part no D1531448



Swivel caster with directional brakes
Helps focus on walking instead of direction. Locks out swivel on caster.
Part no D1531408



Brake system
For the user capable of operating brake function.

Hand-brake kit
Part no D1537656

One-hand brake
Part no Right D1539513
Part no Left D1531450



Oxygen tank holder
Secure mounting of oxygen tank on the walker.
Part no D1531406



Sling
Gives security and makes it possible for the user to start the gait training earlier. No tools needed. Easy to fit and adjust, quick release.
Part no D1536846



Stand Assist includes sling
For additional support while standing. (standing plate + knee support + sling)
Part no D1531449

Standing plate only
Part no D1535330



IV holder
Separate IV-stand makes life complicated. Not anymore.
Part no D1531407



Ball grips
when having difficulties to close the hand.
Part no D1539515

Others

Small caster 4" diameter.
Part no. D153604

Models



Manual
Pull out the pin, move to required position, release the pin.



Pull up Hydraulic
Pull up lever to adjust to required height.

Technical data

Step Up Manual	22.6"/26"/29.6"	27.3"/30.8"/34.3"	34" - 51.5"	13.2" - 18.3"	24" - 29.25"
Step Up Hydraulic	22.6"/26"/29.6"	27.3"/30.8"/34.3"	34" - 51.5"	13.2" - 18.3"	24" - 29.25"

Step Up Manual	31.5"	35.2 lb	330 lb	5"	8.3"
Step Up Hydraulic	31.5"	39.6 lb	330 lb	5"	8.3"

Color



RAL 1035

Product Number

Step Up Manual 1521677
Step Up Hydraulic 1521678



Your US Distributor



Clarke Health Care Products
1003 International Drive
Oakdale PA 15071
www.clarkehealthcare.com
Toll Free: 888-347-4537
Fax: 724-695-2922

Invacare Dolomite AB
Växjövägen 303 75
Box 200, 343 75 Diö
Sweden
Tel. +46 (0)476 535 00
Fax +46 (0)476 535 99
www.dolomite.biz

© 2012 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

StepUp - EU - 09/2012



Yes, you can.®