



Made for training

The base frame is designed to enhance premium gait area. The raised frame in front allows longer step. The width in side the frame generous.

Ergonomic adjustment of arm support

Arm pads are in soft material for pressure relief. Securely assembled with microfilm. Quick and simple to remove for cleaning or replacement.

Transfer and access

Well thought out for transfer situations. From bed, the lowered back leg helps reversing under the bed. When

transferring from wheel chairs benefit of the generous frame size and, the width adjustment.

Features and Options



By releasing the central knob, both arm pads are simultaneously adjusted in width. The position of the arm supports remains at an angle of 25 degrees from the shoulder points - the ideal position for keeping upper strength.

Accessories and options

Quick and easy assembly in track or prepared recess.



Practical basket
Keeping belongings within easy to reach range.
Part no D1531390



Tray
for small things.
Part no D1531391



Elbow stop
Securing elbows from falling out of position.
Part no D1531448



Swivel caster with directional brakes
Helps focusing on walking instead of direction.
Part no D1531408



Brake system
For the user capable of operating brake function.

Hand-brake kit
Part no D1537656

One-hand brake
Part no Right D1539513
Part no Left D1531450



Oxygen tank holder
Secure mounting of oxygen tube on the walker.
Part no D1531406



The sling gives security and make it possible for the user to start the walk training earlier. No tools needed. Easy to fit and adjust, quick release.
Part no D1536846



Stand assist incl. slings
For verticalisation.
Part no D1531449 (standing plate + knee support + sling)
Standing plate
Part no D1535330



IV holder
Separate IV-stand makes life complicated. Not anymore.
Part no D1531407



Ball grips
when having difficulties to close the hand.
Part no D1539515

Others

Small caster 4" diameter.
Part no. D153604

Models



Manual

Pull out the pin move to required position release the pin.



Pull up Hydraulic

Pull up lever to adjust to required height.

Technical data



22.6"/26"/29.6"



27.3"/30.8"/34.3"



34" - 51.5"



13.2" - 18.3"



24" - 29.25"

Step Up Manual

Step Up Hydraulic



31.5"



35.2 lb



330 lb



5"



8.3"

Step Up Manual

Step Up Hydraulic

Colour



RAL 1035

Article code

Step Up Manual 1521677
Step Up Hydraulic 1521678

Invacare Dolomite AB
 Växjövägen 303 75
 Box 200, 343 75 Diö
 Sweden
 Tel. +46 (0)476 535 00
 Fax +46 (0)476 535 99
 www.dolomite.biz

© 2012 Invacare International Sàrl
 All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

StepUp - EU - 09/2012



Yes, you can.®

