

Mobeli Handi-Grip grab bars and accessories

Mobeli is a mobile-holding system that monitors the vacuum state and adhesive strength. Mobeli aids are well-suited for traveling, they can be quickly taken apart, stowed in luggage and easily put in the new bathroom.

The handles have been developed and tested to support natural movement, to stand or sit, or a quick grasp in case of a fall or slip. The capacity indicates the weight tested to break the suction. The red indicators monitor the suction and suggest when to reapply the suction rocker.

These grab handles can be applied to a non-porous surface or tile that is greater than 5"x5". The pad cannot cross grout lines or seams. They are not suitable to mount on new dry wall or wallpaper. Some natural stone tiles may still be porous and thus are not suitable. Allergy notice, the suction pads are made from natural rubber latex.

Product number	Description	Length	Capacity
R1400202S	Shower Head holder, short	9"	
R1499202S	Shower Head holder w/swivel	11.7"	
R1400205S	Shower Screen arm only	28"-10.9"	
R14061208	Curtain only for shower screen	46.8"X31.2"	
R1499205S	Shower Screen and arm combination	11.7"-28.5"	
R1400220S	One hand grip, fixed length	12.8"	231.5 lb
R1400221S	Small telescoping bar	17"-21.4"	187.4 lb
R1400222S	Medium telescoping bar	21.8"-26.7"	187.4 lb
R1400223S	Large telescoping bar	25.7"-30.6"	187.4 lb
R1400224S	X-large telescoping bar	34.5"-39.3"	176.4 lb
R1400241SU	QuattroPlus, 2 swivel points, 4 pads	14"-18.9"	352.8 lb
R1400242SU	QuattroPlus, 2 swivel points, 4 pads	18.7"-24.37"	352.8 lb
R1400243SU	QuattroPlus, 2 swivel points, 4 pads	22.6"-27.49"	352.8 lb
R1400244SU	QuattroPlus, 2 swivel points, 4 pads	21.4"-26.3"	352.8 lb
R1400261SU	QuattroPower Tub, long, 3 swivel points, 4 pads	28.8"-38.2"	176 lb
R1400262SU	QuattroPower Tub, long w/grip horn, 3 swivel points, 4 pads	28.8"-38.2"	176 lb
R1400263S	QuattroPower tub, short, 3 swivel points, 4 pads	28"	176 lb
R1400264S	QuattroPower tub, short w/grip horn, 3 swivel points, 4 pads	33.15"	176 lb
R1400265SU	QuattroPower, 3 swivel points, 4 pads	30.4"	176 lb
R1400266SU	QuattroPower w/grip horn, 3 swivel points, 4 pads	30.4"	176 lb
R1400555	Cardan joint adapt, pr, alum, swing and rotate, decreases capacity on any bar to 110 lb.	3.6"	110 lb
R1400557	Diagonal joint adapt, pr for Quattro bars only, decreases capacity to 165 lb.		165 lb
R1400558	Angle joint adapt, pr, plastic, swings on one plane, reduces capacity to 143 lb		143 lb
R1400581	Dual hook, slides onto any telescoping bar to hold bath tools	3.5" drop	
R1461414	Stainless steel plates, adhesive back, 5.6"x5.6", to secure to surfaces not suitable for suction pads	5.6"x5.56"	

Mobeli Handi-Grip grab bars and accessories

Handi-Grip grab bars with powerful suction pads that allow you to mount a grab bar at the best position for bathing, showering or anywhere at home where you need assistance standing or sitting.




Shower head holders allow a shower wand to be mounted for best position



Shower curtain and mounting arm. For personal privacy.

Grab bars come in fixed lengths and also telescoping lengths for best fit in your shower or tub area.



Quattro bars have 4 suction pads for greater adhesion and capacity. QuattroPower bars have 2 long bars. They can span across a tub or shower corner. Use beside a toilet for assistance.

Joint adapters allow suction pad to swing or swivel to cross corners.

R1400557



R1400555



R1400558

