


Clarke Health Care Products

	Ramp Specifications							
	Personal Mini	Personal	Multipurpose 4'4"	Multipurpose 5'5"	Multipurpose 6'6"	Multipurpose Edgebarrier 4'4"	Multipurpose Edgebarrier 5'5"	Multipurpose Edgebarrier 6'6"
Product number	MP24	MP210	MP44	MP55	MP66	MPEB44	MPEB55	MPEB66
Length	2'4"	2'10"	4'4"	5'5"	6'6"	4'4"	5'5"	6'6"
Width	29"	33"	29.5"	29.5"	29.5"	32"	32"	32"
Folded Dim.	2.75"x11"x28"	2.75"x11"x34"	52"x6"x8"	65"x6"x8"	78"x6"x8"	52"x6.75"x8"	65"x6.75"x8"	78"x6.75"x8"
Weight	8.5 lbs	13 lbs	22 lbs	26.5 lbs	31 lbs	26.5 lbs	31 lbs	37.5 lbs
Capacity	660 lbs	660 lbs	660 lbs	660 lbs	660 lbs	660 lbs	660 lbs	660 lbs
Suggested rise	up to 7"	up to 9"	up to 13"	up to 16"	up to 20"	up to 13"	up to 16"	up to 20"