

Aquatec® | **Bath lift**



Invacare® **H605 Rio**

The *Invacare®* **H605 Rio** offers the ideal support to have a safe bath. With its new design it integrates perfectly in every bathroom. The anodised aluminium frame of the bath lift is easy to clean and to dismantle quickly without any tools.

The bath lift is battery-driven and based on a maintenance-free and robust lead battery.

The stationary backrest offers a comfortable seating position – even without using the seat covers that are additionally available.

AQUATEC®



Yes, you can.®



Invacare® **H605 Rio**

- Stable aluminium frame, easy to clean
- Battery-driven, powerful and silent
- Stationary backrest
- The textured surface and seamless connection between seating surface and side flaps allow use without cover if desired
- Patented, self-release suction cups: can be easily removed
- Removable and machine-washable covers are available as option
- Can be dismantled into four main components

Options and features



Smooth surface: easy to clean.



Easy to disassemble for transport.



Covers available as option.




Height adapters increase the bath lift by 4 cm.

Technical data

Order numbers	Aquatec® H605 Rio
Bath lift without cover	1533957
Cover set	1534619
Height adapters	1534623

Dimensions	
Seat width without side flaps	37,0 cm
Seat width with side flaps (A)	72,0 cm
Seat depth (B)	44,0 cm
Seat height max. (C)	44,0 cm
Seat height min.	6,0 cm
Width of backrest (D)	35,0 cm - 40,0 cm
Height of backrest (E)	61,5 cm
Total depth (F)	68,0 cm
Total height with backrest (G)	103,0 cm
Base area	67,0 cm x 37,0 cm cm
Cover thickness	0,3 cm



More info	
Total weight without battery	9,5 kg
Weight, seat	1,6 kg
Weight, backrest	5,7 kg
Weight battery	1,6 kg
Load capacity	120,0
Colour	White